







CELLO STRIPPER



REPLACEMENT CELLO STRIPS

STRIPS WITH BLANKS	STRIPS WITHOUT BLANKS
 LM 5502 one-sided, 2 mm	 LM 5512 one-sided, 2 mm
 LM 5503 one-sided, 3 mm	 LM 5513 one-sided, 3 mm
 LM 5503D double-sided, 3 mm	 LM 5513D double-sided, 3 mm

The Cello Stripper Holder accepts 2 sizes of steel abrasive strips- 2mm and 3mm; single-sided and double-sided.

KIT CONTAINS:

- Stainless Steel Strip Holder
- Sample Package of 10 Abrasive Strips
(2mm or 3mm width, single-sided or double-sided with and without blanks)