

# DIAMOND SAW & MANDRELS



DIAMOND SAW  
& MANDREL SET

- 2<sup>3</sup>/<sub>4</sub>" diamond coated blade with special safety mandrel
- Cut neat clean dies from casts
- Use on standard dental lathe

**SET CONTAINS:**

Diamond Saw, Mandrel and Wrench  
(All items also available separately)



STAINLESS STEEL  
MANDRELS

Packages of 12 or 100



BRASS HEAD  
MANDRELS

Packages of 12 or 100



NICKEL PLATED  
MANDRELS

Packages of 12 or 100



LARGE HEAD  
MANDRELS

Packages of 12 or 100