

MIXING BOWLS

FLEXIBLE



These flexible, chemical resistant, neoprene rubber mixing bowls are perfect for mixing plaster, dental stone, and silicone casting materials.

After use, just bend the cup and leftover plaster and dental stone simply pops off. Remaining silicone can easily be peeled off after the silicone cures.

AVAILABLE SIZES:

Small- 3¹/₄" W x 2¹/₄" H 150 ml

Medium- 4¹/₄" W x 2³/₄" H 350 ml

Large- 4³/₄" W x 3⁵/₈" H 600 ml

X-Large- 5" W x 4¹/₂" H 850 ml