

# STRIPPING TOOLS

ENAMEL  
HAND STRIPPER



DIAMOND  
DISC GUARD

## ENAMEL HAND STRIPPER KIT

Manual operation means more control and increased safety. The Stainless Steel Enamel Stripper can be used up to the first molar and accepts all sizes of steel abrasive strips and many types of linen, cloth or plastic strips. Single and double sided.

One abrasive strip can be used for several applications. Using the correct width prevents gingival tissue injury.

### KIT INCLUDES:

- 1 Stainless Steel Enamel Hand Stripper
- Abrasive Strips Sample Pack

### ADDITIONAL STRIPS AVAILABLE:

Single sided strips- package of 10 pcs.  
Double sided strips- package of 10 pcs.

### STRIPS AVAILABLE IN:

Fine, Medium and Coarse

### WIDTHS AVAILABLE:

4mm, 6mm, 8mm



## DIAMOND DISC GUARD

Get peace of mind with this invaluable tool by providing the safest possible working environment for your patients. Crafted from stainless steel, the Diamond Disc Guard protects the mouth from spinning Diamond Discs. Fits all 7/8" diamond discs, for use on doriot style handpieces.