

# STRIPPING TOOLS



CONTOURING KIT

## CONTOURING KIT

Removes interproximal tooth mass to produce additional intra-arch space. May also be used for opening contact points, recontouring or removing overhanging fillings. Reusable aluminum-oxide contouring blades.

### KIT INCLUDES:

- 1 Handle
- 12 Assorted Strips

### ADDITIONAL STRIPS AVAILABLE:

Single sided strips- package of 10 pcs.  
Double sided strips- package of 10 pcs.

### STRIPS AVAILABLE IN:

Fine, Medium and Coarse



MINI STRIPPER

## MINI STRIPPER

- Small and easy to handle
- Complete with 16 refill strips (single side only)

### ADDITIONAL STRIPS AVAILABLE:

Serrated strips- package of 10 pcs.  
Steel strips- package of 10 pcs.  
Diamond strips- package of 10 pcs.

