

SLIM-LINE™ VIEW BOX

ULTRA FLAT- ONLY 3/4" THICK!



The attractive Slim Line™ View Box is only 3/4" thick and can easily hang on a wall, lie flat or at an angle with built-in leg support.

Perfect for patient presentations, cephalometric tracing as well as day in and day out use.

Generous 12" x 8" viewing area accommodates a large variety of film sizes- the outer dimensions measure 14 1/2" x 10 3/4".

Cold cathoid technology assures reliability as well as extremely long lasting bulb life.

Unit includes 1 AC Adapter

FULL ONE YEAR WARRANTY