

WAX POTS



ANALOG WAX POT



DIGITAL WAX POT

For the fast and efficient heating of wax. Each Wax Pot features 3 wax wells which allows working on multiple wax types at the same time. Solid, strong construction and 4 non-skid rubber safety stands ensure safety while in use. Their compact size allows each to fit into even the most crowded work spaces.

ANALOG WAX POT

- Temperature dial provides quick control with "Low" and "High" settings
- Temperature Range: 122°F – 239°F
- Efficient: quickly reaches and maintains desired temperature
- Accurate: makes perfect copings every time
- Compartment Depth: 7/8"
- Size: Approx 5" W x 4⁵/₁₆" D x 2¹/₄" H
- Weight: 2lbs
- 110V 60Hz

DIGITAL WAX POT

- Precise digital display of temperature at a glance
- Easy to use "up" and "down" temperature control buttons
- Temperature Range: 122°F – 239°F
- Efficient: quickly reaches & maintains desired temperature
- Accurate: makes perfect copings every time
- Compartment Depth: 5/8"
- Size: Approx. 5¹/₈" W x 6¹/₂" D x 2⁵/₈" H
- Weight: 2.5lbs
- 115V 60Hz